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**Medical device regulation**

This product complies with EU Regulation MDR 2017/745 for medical devices

**1 Purpose**

Electrosurgical equipment for cutting and coagulation of tissue.

**Indications:**

Monopolar forceps are intended to use for electrosurgical coagulation during surgical procedures.

**Contraindications:**

Do not use monopolar forceps in direct contact with the heart, central circulatory system or the central nervous system.

Do not use monopolar forceps where its surgical techniques are contraindicated.

Do not use monopolar forceps if, in the opinion of an experienced physician or according to current professional literature, such use would cause endangerment of the patient, for example due to the general condition of the patient, or if other contraindications are present.

**2 Physical principle of operation**

The flow of current between electrodes / jaws and tissue causes a coagulation effect.

The thermal effect is obtained when you use high-frequency current with a frequency >200 kHz. The electrolytic effect and the irritation effect are not detectable above a frequency of approx. 20 kHz. The coagulation of biological tissue is based on the thermal effect.

**3 Safety instructions**

**3.1 Device-related**

- Goods which are delivered non-sterile must be cleaned and sterilised before use! Sterilisable up to 137°C.
- Do not sterilise the product in hot air.
- Check the device for damage after reprocessing and before use  
**CAUTION:**  
A simple visual inspection on its own is not sufficient to ensure the integrity of the insulation.
- Do not use the product if it is damaged.
- Ensure that the product is cleaned, disinfected and sterilised correctly.
- Do not strain the product beyond its mechanical limits.
- Do not carry out any repair work or maintenance on the device.
- Risk of injury from sharp edges.
- The polished precious metal tips of monopolar forceps may show signs of discolouration. This does not impair their function.
- Please note the instructions for use for the HF generator and connecting cable.
- Sterilise the product before it leaves the hospital/practice environment to prevent the spread of germs and infections, .

**3.2 Use-related**

- Users must be trained in, and be familiar with, the basic principles & rules of the application, and risks of HF surgery.
- The products may only be used at medical facilities by trained medical personnel.
- Ensure that there are no explosive gases or inflammable fluids and materials in the surgical area.
- Report any incidents or near misses to your country's medical devices regulatory agency. When doing so, follow your own facility's internal reporting system. Notify your local BOWA specialist retailer in such cases.
- In the case of patients with pacemakers, consult a cardiologist before the application.
- Comply with the recommended power settings and maximum voltage.
- Always select the minimum required power; the effectiveness of the setting must be assessed by the user.
- Before application, perform a functional test using a piece of gauze soaked with saline solution.
- Coagulate only when the contact surface is within your field of view.
- Not for use for contraceptive tubal coagulation.
- The reliability of coagulation must be assessed by the attending physician according to the nature of the vascularised tissue and the vessel pathology (arteriosclerosis, aneurysms, vessels, etc.).
- Do not coagulate tissue bundles with unknown content.

- Proceed cautiously when performing coagulation in the vicinity of sensitive structures, such as nerves or ureter.
- Avoid contact with metallic objects (clips, stents, etc.) in the area of the active electrode surfaces. This can affect energy output and lead to undesired effects.
- To avoid bodily injuries resulting from inadvertent HF activation, do not place active HF instruments on the patient.
- Make sure that the instrument does not contact other non-insulated instruments or objects (does not apply to BOWA forceps without plug).
- Soiled forceps tips can result in inadequate coagulation.
- The tips may still be hot enough to cause burns immediately after HF power has been switched off.
- To prevent burns from inadvertent activation by foot, the generator must be deactivated or the instrument must be disconnected from the generator for cleaning of the active (connected) forceps.
- The "BF"/"CF" applied part of the HF device being used is extended by any connected instrument.
- In order to prevent electric shock or burns, only contact the forceps within the insulated area.
- To prevent forming a closed circuit when using forceps without a plug, only touch the patient with one hand at a time. Electric shocks may result from a circuit being formed.

### 3.3 EMC information

Medical electrical devices are subject to special EMC precautions, so please follow the instructions detailed below.

- BOWA accessories are only intended for connection to the specified HF devices.
- Using the accessory with medical devices from other manufacturers may result in higher emission levels or reduced interference immunity.

Combining medical devices is safe only if <ul style="list-style-type: none"> <li>➤ the relevant instructions for use permit the desired combination or</li> <li>➤ the intended use and the interface specification of the products used in the combination are compatible.</li> </ul> The instructions for use and interface specifications of medical devices used in combination must be strictly observed.

## 4 Procedure during use

1. Sterilise the forceps, which are delivered in non-sterile condition, before use.
2. Connect the correctly sized forceps to a connection cable that is compatible with the HF generator (see section 4.2 *Connection combinations*). The use of BOWA connection cables is recommended, since the plug connections of other manufacturers might not comply with the technical standards.
3. Do not exceed the maximum permissible voltage.

### 4.1 HF unit/power settings

The maximum permissible voltage is 3500 Vp.  
The instrument must be used with monopolar coagulation current exclusively.

### 4.2 Connection combinations

For permissible connection combinations, see the BOWA compatibility declaration.

## 5 Reprocessing

### 5.1 General information

The following values for the permitted number of reprocessing cycles are intended as a guide. The actual number may vary depending on the demands placed on the devices.  
Manual reprocessing methods are not recommended due to their significantly lower effectiveness.  
Follow the chemical manufacturers' statements on concentrations and exposure times. The medical device manufacturer cannot be held liable if other types of cleaning agents and disinfectants are used.

### 5.2 Product-specific information

Insulated parts must not come into contact with hard, sharp or heavy products, as this could damage the electrical insulation.  
Do not use wire brushes, abrasives or other products which could damage the surface.

<b>Reprocessing cycles</b>	75
<b>Reprocessing step</b>	<b>Description</b>
<b>Manual pre-cleaning</b>	
<b>1.A Soaking<sup>1</sup></b>	Soak the product at room temperature (<25°C) for at least 5 minutes immediately, but no later than 2 hours, after use. Only use enzymatic, aldehyde-free disinfectants which are suitable for disinfecting medical devices (e.g., DGHM or FDA approval or CE marking). Remove all visible soiling with a soft plastic brush. BOWA recommends the use of Gigasept® Instru AF (Schülke & Mayr GmbH).
<b>1.B Ultrasound<sup>2</sup></b>	Place the device in the ultrasonic bath for at least 5 minutes. BOWA recommends the use of Gigasept® Instru AF (Schülke & Mayr GmbH).
<b>2 Rinsing</b>	Rinse the product thoroughly for at least 1 minute at room temperature (<25°C) under running tap water (at least drinking water quality). Allow residual water to drip off sufficiently.
<b>Automatic cleaning and disinfection</b>	
<b>3 Automatic cleaning<sup>3</sup></b>	Use a washer-disinfector (WD) with tested effectiveness (according to ISO 15883) and a neutral to alkaline (max. pH 11.5) enzymatic cleaning agent. Depending on the concentration, agents including alcohol- and/or aldehyde-containing ingredients may be used. BOWA recommends the use of neodisher® MediClean forte (Chemische Fabrik Dr. Weigert) at 55°C for 10 minutes. The intermediate rinsing must include at least two rinsing steps with fully demineralised water (deionised water) (at least 1 minute each at >10°C).
<b>4 Thermal disinfection<sup>4</sup></b>	An A <sub>0</sub> value of >3000 must be attained. BOWA recommends a temperature of 90°C for at least 5 minutes and the use of demineralised water.
<b>5 Drying</b>	Drying is carried out according to the WD programme and depends on the total load. The maximum temperature is 100°C for 25 minutes. If necessary, dry with filtered compressed air at max. 3 bar.
<b>Inspection</b>	After machine cleaning and disinfection, carry out a visual inspection to check for residues. Repeat the reprocessing process if necessary.
<b>Steam sterilisation</b>	
<b>6 Steam sterilisation<sup>5</sup></b>	BOWA recommends steam sterilisation using the fractionated vacuum method with a sterilisation time of 3–20 minutes and a temperature of 134–137°C in a suitable sterilisation packaging.

The manufacturer does not accept any responsibility if other types of cleaning agents and disinfectants are used. Comply with the recommendations of the manufacturers of the cleaning agents. Please note also the provisions of the sterilisation service provider.

#### The following have been validated:

- 1 Soaking for 5 minutes with Gigasept® Instru AF (Schülke & Mayr GmbH), 3% v/v, brushing
- 2 Ultrasonic bath for 5 minutes with Gigasept® Instru AF (Schülke & Mayr GmbH), 3% v/v
- 3 Neodisher® MediClean forte (Chemische Fabrik Dr. Weigert), 0.5% v/v
- 4 Thermal disinfection at 90°C for 5 minutes
- 5 Half-cycle with 3 pre-vacuum phases at 132°C for 1.5 minutes, packaging: double-packaged in paper/film sterilisation pouch

## 6 Storage and transport

Store the devices in their original packaging and protect them against direct sunlight and humidity.

## 7 Disposal

Dispose of the medical devices, packaging material and accessories in accordance with the respectively applicable national regulations and laws.

## 8 Symbols to bear in mind

Symbol	Meaning	Symbol	Meaning
	Follow the instructions for use		Keep away from direct sunlight
<b>REF</b>	Catalogue number		Keep dry
<b>LOT</b>	Batch code		Conforms to Russian standards
<b>QTY</b>	Quantity		Manufacturer
<b>CE</b>	CE mark	<b>CE</b> 0123	CE mark and identification number of Notified Body (class IIb)
<b>R<sub>Only</sub></b>	Caution: Federal law restricts this device to sale by or on the order of a physician. For the attending physician only.		WEEE symbol: Requires separate collection for waste of electric and electronic equipment.
<b>MD</b>	Medical device		Caution

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